

| SystemRoMedic™

# Lifting Slings



# HighBackSling

## Standard model for users who need a lot of support.

HighBackSling is designed for use in most common lifting situations and when the user needs extra head support. HighBackSling can also be used for lifting from a recumbent position on the floor.

- Divided leg supports; easy for the caregiver to apply and remove, the user does not sit on the lifting sling
- High back; the user is supported from the tailbone up over the head and under the legs
- Slightly reclined seated posture
- Standard model, suits most users
- Removable plastic stays give the head support extra stability
- HighBackSling in polyester net is available with or without reinforced leg supports
- Pocket at lower edge for easier application
- SWL 300 kg/660 lbs (disposable material: SWL 200 kg/440 lbs)

**Accessories:** AdjustmentLoop, SoftLegSupport, StrapsPadding, NameTag

Material / Size



# Size guide for lifting slings

## Divided leg supports / High back

Product name	Size	Rec. user weight	Sling seat measurement	Sling overall height
HighBackSling	XXS	10-25 kg / 22-55 lbs	20 cm / 8"	77 cm / 30"
HighBackSling	XS	20-35 kg / 44-77 lbs	25 cm / 10"	85 cm / 33"
HighBackSling	S	30-55 kg / 66-122 lbs	35 cm / 14"	91 cm / 36"
HighBackSling	M	50-90 kg / 110-198 lbs	45 cm / 18"	94 cm / 37"
HighBackSling	L	80-110 kg / 176-243 lbs	56 cm / 18"	98 cm / 39"
HighBackSling	XL	110-160 kg / 243-353 lbs	62 cm / 24"	98 cm / 39"
HighBackSling	XXL	130-200 kg / 287-441 lbs	70 cm / 28"	102 cm / 40"

ShellHBSling	XXS	10-30 kg / 22-66 lbs	17 cm / 7"	80 cm / 31"
ShellHBSling	XS	20-35 kg / 44-77 lbs	25 cm / 10"	85 cm / 33"
ShellHBSling	S	30-55 kg / 66-121 lbs	35 cm / 14"	94 cm / 37"
ShellHBSling	M	50-90 kg / 110-198 lbs	40 cm / 16"	97 cm / 38"
ShellHBSling	L	80-120 kg / 176-264 lbs	45 cm / 18"	98 cm / 39"
ShellHBSling	XL	110-160 kg / 243-353 lbs	55 cm / 22"	100 cm / 39"
ShellHBSling	XXL	150-210 kg / 331-463 lbs	65 cm / 26"	100 cm / 39"

FlexibleSling	XXS	10-30 kg / 22-66 lbs	20 cm / 8"	75 cm / 30"
FlexibleSling	XS	20-35 kg / 44-77 lbs	25 cm / 10"	81 cm / 32"
FlexibleSling	S	25-55 kg / 55-121 lbs	35 cm / 14"	90 cm / 35"
FlexibleSling	M	50-75 kg / 110-165 lbs	45 cm / 18"	102 cm / 40"
FlexibleSling	L	70-110 kg / 154-242 lbs	55 cm / 22"	102 cm / 40"
FlexibleSling	XL	100-130 kg / 220-287 lbs	65 cm / 26"	104 cm / 41"

## Divided leg supports / Shoulder-high back

Product name	Size	Rec. user weight	Sling seat measurement	Sling overall height
BasicSling	XXS	10-25 kg / 22-55 lbs	21 cm / 8"	63 cm / 25"
BasicSling	XS	20-35 kg / 44-77 lbs	27 cm / 11"	65 cm / 26"
BasicSling	S	35-55 kg / 77-121 lbs	37 cm / 15"	77 cm / 30"
BasicSling	M	45-80 kg / 99-176 lbs	50 cm / 20"	89 cm / 35"
BasicSling	L	75-110 kg / 165-243 lbs	60 cm / 24"	95 cm / 37"
BasicSling	XL	100-150 kg / 220-301 lbs	65 cm / 26"	98 cm / 39"
BasicSling	XXL	140-200 kg / 309-441 lbs	75 cm / 30"	98 cm / 39"

AmpSling	XS	20-40 kg / 44-88 lbs	27 cm / 11"	85 cm / 33"
AmpSling	S	35-60 kg / 77-133 lbs	37 cm / 15"	90 cm / 35"
AmpSling	M	50-80 kg / 110-176 lbs	45 cm / 18"	97 cm / 38"
AmpSling	L	70-110 kg / 154-242 lbs	50 cm / 20"	100 cm / 39"
AmpSling	XL	100-130 kg / 220-287 lbs	60 cm / 24"	100 cm / 39"
AmpSling	XXL	120-180 kg / 265-397 lbs	70 cm / 28"	100 cm / 39"

## Divided leg supports / Low back

Product name	Size	Rec. user weight	Sling seat measurement	Sling overall height
LowBackSling	XXS	10-25 kg / 22-55 lbs	21 cm / 8"	60 cm / 24"
LowBackSling	XS	20-35 kg / 44-77 lbs	27 cm / 11"	60 cm / 24"
LowBackSling	S	30-50 kg / 66-110 lbs	37 cm / 15"	77 cm / 30"
LowBackSling	M	45-80 kg / 99-176 lbs	47 cm / 18"	85 cm / 33"
LowBackSling	L	75-110 kg / 165-243 lbs	57 cm / 22"	90 cm / 35"
LowBackSling	XL	100-140 kg / 220-309 lbs	67 cm / 26"	92 cm / 36"
LowBackSling	XXL	130-190 kg / 287-419 lbs	77 cm / 30"	92 cm / 36"

## Undivided leg supports / High and shoulder-high back

Product name	Size	Rec. user weight	Sling seat measurement	Sling overall height
FlexibleSling, undivided	XS	20-40 kg / 44-88 lbs	75 cm / 30"	80 cm / 31"
FlexibleSling, undivided	S	30-70 kg / 66-154 lbs	85 cm / 33"	85 cm / 33"
FlexibleSling, undivided	M	60-90 kg / 132-198 lbs	94 cm / 37"	93 cm / 37"
FlexibleSling, undivided	L	80-130 kg / 176-287 lbs	105 cm / 41"	96 cm / 38"
FlexibleSling, undivided	XL	110-180 kg / 242-397 lbs	115 cm / 45"	98 cm / 39"

ClassicHBSling	XS	20-40 kg / 44-88 lbs	30 cm / 12"	95 cm / 37"
ClassicHBSling	S	30-60 kg / 66-132 lbs	40 cm / 16"	98 cm / 39"
ClassicHBSling	M	55-85 kg / 121-187 lbs	48 cm / 19"	106 cm / 42"
ClassicHBSling	L	80-120 kg / 176-265 lbs	65 cm / 26"	110 cm / 43"
ClassicHBSling	XL	110-140 kg / 242-309 lbs	80 cm / 31"	115 cm / 45"

ClassicSling	XS	20-40 kg / 44-88 lbs	30 cm / 12"	80 cm / 31"
ClassicSling	S	30-60 kg / 66-132 lbs	38 cm / 15"	90 cm / 35"
ClassicSling	M	55-85 kg / 121-187 lbs	46 cm / 18"	95 cm / 37"
ClassicSling	L	80-120 kg / 176-265 lbs	63 cm / 25"	98 cm / 39"
ClassicSling	XL	110-140 kg / 242-309 lbs	80 cm / 31"	100 cm / 39"

BariSling	XXXL	200-570 kg / 441-1257 lbs	-	-
-----------	------	---------------------------	---	---

## For toileting / High back

Product name	Size	Rec. user weight	Sling seat measurement	Sling overall height
HygieneHBSling	XXS	10-30 kg / 22-66 lbs	25 cm / 9.8"	92 cm / 36"
HygieneHBSling	XS	25-45 kg / 55-99 lbs	28 cm / 11"	92 cm / 36"
HygieneHBSling	S	40-60 kg / 88-132 lbs	30 cm / 12"	96 cm / 38"
HygieneHBSling	M	55-80 kg / 121-176 lbs	35 cm / 14"	102 cm / 40"
HygieneHBSling	L	75-110 kg / 165-220 lbs	41 cm / 16"	102 cm / 40"
HygieneHBSling	XL	90-130 kg / 198-287 lbs	46 cm / 18"	104 cm / 41"
HygieneHBSling	XXL	120-170 kg / 265-375 lbs	60 cm / 24"	104 cm / 41"

FlexibleSling, undivided*	XS	20-40 kg / 44-88 lbs	75 cm / 30"	80 cm / 31"
FlexibleSling, undivided*	S	30-70 kg / 66-154 lbs	85 cm / 33"	85 cm / 33"
FlexibleSling, undivided*	M	60-90 kg / 132-198 lbs	94 cm / 37"	93 cm / 37"
FlexibleSling, undivided*	L	80-130 kg / 176-287 lbs	105 cm / 41"	96 cm / 38"
FlexibleSling, undivided*	XL	110-180 kg / 242-397 lbs	115 cm / 45"	98 cm / 39"

\*) with hygiene opening