

| SystemRoMedic™

# Lifting Slings

# StretcherSling



## For horizontal lifting with StretcherBar.

StretcherSling is a lifting sling that is used in combination with mobile lifts or stationary lift systems and SystemRoMedic's StretcherBar, a specially designed sling bar for horizontal lifting and transferring users in a lying position. StretcherSling has ten suspension points; five on either side of the user, giving the user evenly distributed support under the entire body.

- Easy for the caregiver to apply and remove
- Gives the user an added sense of security
- Five lift straps on both sides of the lifting sling
- StretcherSling is available in two different materials and in one size
- SWL 300 kg/660 lbs

**Accessories:** AdjustmentLoop, NameTag

Material / Size

 **Polyester**  
One size,  
SWL 300 kg/660 lbs

 **Disposable material**  
One size,  
SWL 300 kg/660 lbs



# SafeHandlingSheet



For horizontal lifting and positioning with a two-point sling bar.

SafeHandlingSheet is a sheet that is used for horizontal lifting or for positioning a user in bed without subjecting the caregiver to unnecessary strain. SafeHandlingSheet is used in combination with mobile lifts or stationary lift systems and a two-point sling bar. SafeHandlingSheet has sixteen lift straps, eight on either side of the user, with three alternative lift loops each, which enables optimal adaptation to the user and transfer situation. The lift loops are connected to the sling bar and when the user has been positioned, the lift straps are simply tucked under the mattress or placed in the pockets at the sides of the sheet.

- Eight lift straps, each with three alternative lift loops on either side of the sheet
- Lift loops are colour-coded for safe and easy application
- Pockets for practical storage of lift straps
- Positioning handles at the sides for easy, manual positioning
- SafeHandlingSheet is available in five different materials and sizes to suit different patients and situations

**Accessories:** AdjustmentLoop (SWL 300 kg/660 lbs), NameTag

Material / Size

**Polyester net**  
SWL 300 kg/660 lbs

**Polyester net/WendyLett4Way\***  
SWL 300 kg/660 lbs

**Polyester**  
SWL 400 kg/900 lbs

**Polyester, reinforced**  
SWL 570 kg/1250 lbs



**Disposable material**  
SWL 300 kg/660 lbs



\*) WendyLett4Way is SystemRoMedic's unique sliding drawsheet. Here, we have used the same comfortable, durable and functional material to allow SafeHandlingSheet to slide freely in four directions, without friction or shearing and with minimal effort. Read more about our WendyLett series at [www.handicare.com](http://www.handicare.com).

### For toileting / Low back



Product name	Size	Rec. user weight	Sling seat measurement	Sling overall height
HygieneSling	XXS	10-30 kg / 22-66 lbs	25 cm / 10"	50 cm / 20"
HygieneSling	XS	25-45 kg / 55-99 lbs	35 cm / 14"	65 cm / 26"
HygieneSling	S	40-60 kg / 88-132 lbs	48 cm / 19"	75 cm / 30"
HygieneSling	M	55-85 kg / 121-187 lbs	55 cm / 22"	87 cm / 35"
HygieneSling	L	75-100 kg / 165-220 lbs	60 cm / 24"	93 cm / 37"
HygieneSling	XL	90-130 kg / 198-287 lbs	80 cm / 31"	95 cm / 37"
HygieneSling	XXL	120-170 kg / 265-376 lbs	90 cm / 35"	95 cm / 37"



HygieneLBSling	XXS	10-30 kg / 22-66 lbs	25 cm / 10"	92 cm / 36"
HygieneLBSling	XS	25-45 kg / 55-99 lbs	28 cm / 11"	92 cm / 36"
HygieneLBSling	S	40-60 kg / 88-132 lbs	30 cm / 12"	96 cm / 38"
HygieneLBSling	M	55-80 kg / 121-176 lbs	35 cm / 14"	102 cm / 40"
HygieneLBSling	L	75-110 kg / 165-220 lbs	41 cm / 16"	102 cm / 40"
HygieneLBSling	XL	90-130 kg / 198-287 lbs	46 cm / 18"	104 cm / 41"
HygieneLBSling	XXL	120-170 kg / 265-376 lbs	60 cm / 24"	104 cm / 41"

### Horizontal lifting / Turning



Product name	Size	Rec. user weight	Length	Width
StretcherSling	one size	50-200 kg / 110-441 lbs	180 cm / 71"	86 cm / 34"
StretcherSling, disposable material	one size	50-200 kg / 110-441 lbs	180 cm / 71"	86 cm / 34"



SafeHandlingSheet	L	50-300 kg / 110-661 lbs	240 cm / 94"	115 cm / 45"
SafeHandlingSheet	XL	50-400 kg / 110-881 lbs	240 cm / 94"	115 cm / 45"
SafeHandlingSheet	XXL	50-570 kg / 110-1257 lbs	240 cm / 94"	130 cm / 51"
SafeHandlingSheet, disposable material	M	50-300 kg / 110-661 lbs	200 cm / 79"	120 cm / 47"

### Rising/ Standing with MiniLift



Product name	Size	Rec. user weight	Lift vest width
ThoraxSling	XS	20-50 kg / 44-110 lbs	85 cm / 33"
ThoraxSling	S	30-60 kg / 66-132 lbs	94 cm / 37"
ThoraxSling	M	50-80 kg / 110-176 lbs	104 cm / 41"
ThoraxSling	L	70-110 kg / 154-220 lbs	124 cm / 49"
ThoraxSling	XL	100-150 kg / 220-330 lbs	146 cm / 57"
ThoraxSling	XXL	140-200 kg / 309-441 lbs	170 cm / 67"

### Rising/ Standing with RiseUp600EE



Product name	Size	Rec. user weight	Lift vest width
StandUpSling	S	40-65 kg / 88-143 lbs	109 cm / 43"
StandUpSling	M	55-80 kg / 121-176 lbs	121 cm / 48"
StandUpSling	L	75-110 kg / 165-220 lbs	122 cm / 48"
StandUpSling	XL	100-160 kg / 220-353 lbs	136 cm / 54"
StandUpSling	XXL	150-220 kg / 331-485 lbs	166 cm / 65"

### Standing / Walking



Product name	Size	Rec. user weight	Lift vest width	Lift vest overall height
StandingVest	XS	15-40 kg / 33-88 lbs	82 cm / 32"	68 cm / 27"
StandingVest	S	20-50 kg / 44-110 lbs	96 cm / 38"	73 cm / 29"
StandingVest	M	40-60 kg / 88-132 lbs	110 cm / 43"	73 cm / 29"
StandingVest	L	50-80 kg / 110-176 lbs	122 cm / 48"	73 cm / 29"
StandingVest	XL	70-100 kg / 154-220 lbs	130 cm / 51"	73 cm / 29"



WalkingVest	XXS	15-40 kg / 33-88 lbs	60 cm / 24"	67 cm / 26"
WalkingVest	XS	25-50 kg / 55-110 lbs	70 cm / 28"	68 cm / 27"
WalkingVest	S	40-60 kg / 88-132 lbs	82 cm / 32"	68 cm / 27"
WalkingVest	M	55-80 kg / 121-176 lbs	96 cm / 37"	73 cm / 29"
WalkingVest	L	75-100 kg / 165-220 lbs	110 cm / 43"	73 cm / 29"
WalkingVest	XL	90-110 kg / 198-243 lbs	122 cm / 48"	73 cm / 29"
WalkingVest	XXL	100-130 kg / 220-287 lbs	142 cm / 56"	75 cm / 30"



BariVest	XXXL	180-300 kg / 397-661 lbs	190-280 / 75-110"	
BariVest	XXXXL	250-570 kg / 551-1257 lbs	230-310 / 91-122"	

### Please note!

To determine the most appropriate model, size and material for the individual user, trial fitting of lifting slings is always essential.

The size guide's measurements and recommendations are only intended as an aid to selecting lifting slings.

The lifting sling's seat width gives an indication of how wide the user's hips (trochanter width) can be. The user's trochanter width must be within +/- 10 cm/3.9" of the lifting sling's seat width.

The lift vest's width gives an indication of the girth of the user's trunk. The user's trunk girth must be within +/- 10 cm/3.9" of the lift vest's width.

The lifting sling is too large if the user's seat slides out too far at the lower opening, or if the leg supports slide forward towards the hollows of the user's knees.

The lifting sling is too small if the user's face is too close to the slingbar, if the leg supports chafe at the groin, or if the sling appears to be too tight or feels too tight.

